DNA ekstraksjon gruppe 2

Date: 2024-09-05

Tags: DNA

Created by: Sivert Helland Veseth

Purpose:

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Extraction of DNA from whole blood / Utført 05.09.2024

Background

This protocol is adopted from Bartlett, J. M. S. and D. Stirling (2003), PCR protocols, Ch 6.

Materials

- 15 ml centrifuge tubes
- 1.5 ml centrifuge tubes
- Waterbath (65°C)
- Tube rotator

Methods

Blood can be collected in heparin or EDTA-vacutainer tubes. Process samples on the same days as collection if stored at room temp.

1. Move 3mL of whole blood to a 15 mL tube.
2. Add 12 mL of **Reagent A**
3. Mix slowly by rotation, 4 min at room temp.
4. Centrifuge 3000 g for 5 min at room temp.
5. Discard supernatant without disturbing cell pellet. Remove all moisture by inverting and blotting on paper or by pipetting.
6. Add 1 mL of **Reagent B**, vortex to re-suspend the cell pellet.
7. Add 250 µl of **5 M sodium perchlorate** and mix by inverting tube several times. ⚠ Sodium perchlorate is harmful if swallowed and causes irritation. Keep away from heat.
8. Place the tube in 65°C waterbath for 15-20 min.
9. Cool to room temp. 10 Add 2 mL of ice-cold **chloroform**. ⚠ Toxic, perform operation in fume hood with appropriate protective equipment.
10. Mix on a rotating mixer for 30 to 60 minutes (more or less will reduce DNA yield).
11. Centrifuge, 2400g for 2 min. Careful not to disturb the separated phases.

12. Transfer upper phase into a clean falcon tube with a sterile pipette.
13. Precipitate the DNA by adding 2-3 mL of chilled **100% ethanol**, invert tube. DNA should appear in the solution.
14. Transfer the DNA to a 1.5 mL tube allow to air dry.
15. Resuspend in 200 ul of **TE buffer**.
16. Quantify DNA concentration with the spectrophotometer, 200-500 ng/ul is expected.

Solutions

0.5 M EDTA, pH 8, for 1 l:

COMPONENT	WEIGHT
EDTA	146.1 g
NaOH	~20g

Add EDTA to 800 mL ddH₂O, adjust pH with NaOH. Autoclave 15 psi, 15 min.

1 M Tris-HCl, pH 7.6, for 1 L:

COMPONENT	WEIGHT
Tris-base	121.1 g

Add to 800 mL ddH₂O, adjust pH with HCl and make up to 1L. Autoclave 15 psi, 15 min.

Reagent A, red blood cell lysis, pH 8. For 1L:

COMPONENT	CONCENTRATION	WEIGHT/VOLUME
Tris-HCl	0.01 M	10 ml 1 M Tris-HCl
Sucrose	320 mM	109.54 g
MgCl ₂	5 mM	0.47 g
Triton X-100	1%	10 mL

Adjust pH to 8 with NaOH. Autoclave 10 psi, 10 min.

Reagent B, cell lysis

COMPONENT	CONCENTRATION	WEIGHT/VOLUME
Tris-HCl	0.4 M	400 mL 1 M Tris-HCl
NaCl	150 mM	8.76 g
EDTA	0.06 M	120 mL 0.6 M EDTA
SDS	1%	10 g ¹

Adjust pH to 8 and make to 1 L with ddH₂O. Autoclave at 15 min at 15 psi².

5 M sodium perchlorate, for 100 mL:

COMPONENT	WEIGHT/VOLUME
sodium perchlorate	70 g

Dissolve in 80 mL ddH₂O make to 100 mL.

TE-buffer, for 1 L:

COMPONENT	CONCENTRATION	WEIGHT/VOLUME
Tris-HCl	0.01 M	10 mL 1 M Tris-HCl ph 7.6
EDTA	1 mM	2 mL 0.5 EDTA

Adjust pH to 7.6, make to 1 L, autoclave 15 min at 15 psi.

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Determination of ACTN3 genotype / Utført 06.09.2024

Background

This protocols is adopted from Schadock et al. (2015)¹ and is suitable for determining the ACTN3 genotype.

Materials

- A PCR machine
- Specific primers

Primers used un Schadock et al.

PRIMER	SEQUENCE	FINAL REACTION CONCENTRATIONS	VOLUME OF 100 UM IN 500 UL PRIMER MIX
hsACTN3_F1	CGCCCTTCAACAACTGGCTGGA	0.5 uM	10 ul
hsACTN3_R1	GATGAGCCCGAGACAGGCAAGG	0.5 uM	10 ul
hsACTN3Tif_F2	CAACACTGCCCGAGGCTGACTG	0.125 uM	2.5 ul
hsACTN3Cir_R2	CATGATGGCACCTCGCTCTCGG	0.25 uM	5 ul

Add primers to 472.5 ul nuclease free H₂O to make a 4X primer mix.

- PCR master mix
- PCR reactions tubes

Methods

- Extract DNA (from whole blood or muscle)
- Mix sample, master mix and primer mix per reaction in a PCR reaction tube:

COMPONENT	VOLUME
2X Master mix	10 ul
Primer mix	5 ul
Sample (correspondning to 100-250 ng genomic DNA)	5 ul

- Run the PCR reaction

TIME	TEMPERATURE	CYCLES
2 min	95°C	
10 sec	95°C	
10 sec	68°C	35 cycles
45 sec	72°C	
2 min	72°C	

- Run the PCR product on a 2% agarose gel using 2-5 ul of the reaction mix. Use a DNA ladder and H₂O as negative control.

Analysis

A R/R genotype produces two bands at 690 and 413 bp. A X/X genotype produces two bands at 690 and 318 bp.

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General DNA Agarose Gel Electrophoresis protocol

Background

This protocol is suitable for electrophoresis separation of DNA fragments from PCR reactions.

Solutions

- 10X TBE, for 1 l:

COMPONENT	WEIGHT	CONCENTRATION
Tris base	121.1 g	1M
Boric acid	61.8 g	1 M
EDTA (disodium salt)	7.4 g	0.02 M

- 6X DNA loading dye, for 10 ml:

COMPONENT	AMOUNT	CONCENTRATION
Glycerol	3 ml	30 %
Bromophenol blue	25 mg	0.25%

Make to 10 ml final volume.

Materials

- Horizontal electrophoresis unit with power supply
- Conical flask for gel preparations
- Heat (microwave or hotplate)
- Gel combs

Methods

Gel preparation

- Dilute 10X TBE buffer to make a 1X solution (100 ml in 900 ml H₂O).
- Add 100 ml of 1X TBE to a conical beaker. Add agarose to create a suitable percentage gel (2 g for 2% gel).
- If using Sybr Safe gel stain, add the stain 1:10000 to the mixture (10 ul in 100 ml)
- Heat the mixture until the solution gets clear (e.g. 1 min in microwave or on a hotplate).
- Let the mixture cool to about 60 degrees (it should feel warm but not hot).
- Pour the mixture into a gel casting tray and place the comb at its designated place.
- The gel will polymerize in 1 h, remove the comb after polymerization.

Running the gel

- Place the gel in the electrophoresis unit.
- Pour 1X TBE into the electrophoresis tank to cover the well.
- Mix samples with loading dye (1 ul of 6X dye per every 5 ul of sample).
- Add DNA ladder and samples to the wells
- Set voltage to 150 V, run four about 1 h (until dye is about 80% through), check that samples moves towards the positive elctrode!

Visualisation

- Visuialize the gel i G:Box using the UV light emitter and Sybr green settings.

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*Wilfinger WW, Mackey K, Chomczynski P. Effect of pH and Ionic Strength on the Spectrophotometric Assessment of Nucleic Acid Purity. BioTechniques 1997;22(3):474-81.

*<https://assets.thermofisher.com/TFS-Assets/CAD/Product-Bulletins/TN52646-E-0215M-NucleicAcid.pdf>

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Attached files

bloxberg-proof_2024-09-06T1318420200.zip (Timestamp archive by Sivert Helland Veseth)
sha256: 2dabf449b84660ee44c07653010c0690b74fab58800931ad32863d09a2cc368a

IDR4000_DNA_agarose_ACTN3_gr.12.tif (Gr. 2 har brønn 2 til 8.)
sha256: a1572bc77bfc150c984ef92d37ed8d801bcd67b27a8e82c9cecc093b477829b

Unique eLabID: 20240905-10827dc6fdeb94d74d5c478021ab156466c35cb7
Link: <https://elab.inn.no/experiments.php?mode=view&id=84>